

Product Name: Septin 8 Rabbit Monoclonal antibody
Catalog #: AMRe03199



Summary

Production Name	Septin 8 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Hamster,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	SEPTIN8
Alternative Names	SEP2; SEPT8
Gene ID	23176
SwissProt ID	Q92599.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 56 kDa; Observed MW: 56 kDa

Background

Filament-forming cytoskeletal GTPase (By similarity).


Product Name: Septin 8 Rabbit Monoclonal antibody
Catalog #: AMRe03199



Research Area

-

Image Data



Western blot analysis of Septin 8 in CHO-K1, C6 lysates using Septin 8 antibody.

Note

For research use only.