

Product Name: SMARCD1 Rabbit Monoclonal Antibody**Catalog #: AMRe03207**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 58 kDa; Observed MW: 58 kDa

Antigen Information

Gene Name	SMARCD1
Alternative Names	CSS11; Rsc6p; BAF60A; CRACD1
Gene ID	6602
SwissProt ID	Q96GM5
Immunogen	A synthetic peptide of human SMARCD1


Background

Involved in transcriptional activation and repression of select genes by chromatin remodeling (alteration of DNA-nucleosome topology).

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of SMARCD1 in HeLa, A549, HL-60 lysates using SMARCD1 antibody.