

Summary

Production Name	STAT6 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	STAT6
Alternative Names	STAT6; Signal transducer and activator of transcription 6; IL-4 Stat
Gene ID	6778
SwissProt ID	P42226.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 94 kDa; Observed MW: 110 kDa

Background

STAT6 transcription factor of the STAT family. Plays a central role in IL4-mediated biological responses. Induces the

Product Name: STAT6 Rabbit Monoclonal Antibody
Catalog #: AMRe03215




expression of BCL2L1/BCL-X(L), which is responsible for the anti-apoptotic activity of IL4.

Research Area

Signal Transduction

Image Data



Western blot analysis of STAT6 in HeLa, U2OS lysates using STAT6 antibody.

Note

For research use only.