

Product Name: Transcription Initiation Factor IIB Rabbit Monoclonal Antibody**Catalog #: AMRe03228**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.68mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 35 kDa; Observed MW: 35 kDa

Antigen Information

Gene Name	GTF2B General transcription factor IIB; General transcription factor TFIIIB; gtf2b; RNA polymerase II
Alternative Names	transcription factor IIB; S300 II; S300-II; TF IIB; TF2B; TF2B_HUMAN; TFIIIB; Transcription initiation factor IIB.
Gene ID	2959
SwissProt ID	Q00403
Immunogen	A synthetic peptide of human TFIIIB


Background

General factor that plays a major role in the activation of eukaryotic genes transcribed by RNA polymerase II.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of TFIIIB in 3T3, Raw264.7 lysates using Transcription Initiation Factor IIB antibody.



Western blot analysis of TFIIIB in HeLa, A549, HL-60, U2OS, C6 lysates using TFIIIB antibody.