

Summary

Production Name	WWOX Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	WWOX
Alternative Names	FOR; WOX1; EIEE28; FRA16D; SCAR12; HHCMA56; PRO0128; SDR41C1; D16S432E
Gene ID	51741
SwissProt ID	Q9NZC7.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 47 kDa; Observed MW: 47 kDa

Background

Product Name: WWOX Rabbit Monoclonal Antibody
Catalog #: AMRe03256




Putative oxidoreductase.


Research Area

Cell Biology


Image Data



Western blot analysis of WWOX in mouse testis, mouse spleen lysates using WWOX antibody.



Western blot analysis of WWOX in HeLa, U2OS, U251MG lysates using WWOX antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using WWOX antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Product Name: WWOX Rabbit Monoclonal Antibody
Catalog #: AMRe03256



Note

For research use only.