
Product Name: DAP kinase 3 Rabbit Monoclonal Antibody**Catalog #: AMRe03262**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 53 kDa; Observed MW: 53 kDa

Antigen Information

Gene Name	DAPK3
Alternative Names	DAPK3; ZIPK; Death-associated protein kinase 3; DAP kinase 3; DAP-like kinase; Dlk; MYPT1 kinase; Zipper-interacting protein kinase; ZIP-kinase
Gene ID	1613
SwissProt ID	O43293
Immunogen	A synthetic peptide of human ZIP Kinase

Background


Death-associated protein kinase 3 (DAPK3) induces morphological changes in apoptosis when overexpressed in mammalian

cells. These results suggest that DAPK3 may play a role in the induction of apoptosis.

Research Area

Cell Biology

Image Data



Western blot analysis of ZIP Kinase in HeLa, A549 lysates using DAP kinase 3 antibody.