

Product Name: ZNF217 Rabbit Monoclonal antibody**Catalog #: AMRe03263**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Concentration	0.47mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 115 kDa; Observed MW: 115 kDa

Antigen Information

Gene Name	ZNF217
Alternative Names	ZABC1; ZN217; ZNF217
Gene ID	7764
SwissProt ID	O75362
Immunogen	A synthetic peptide of human ZNF217


Background

Binds to the promoters of target genes and functions as repressor. Promotes cell proliferation and antagonizes cell death. Promotes phosphorylation of AKT1 at 'Ser-473'.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of ZNF217 in A549, HL-60, U251 lysates using ZNF217 antibody.