
Product Name: MHC Class I Rabbit Monoclonal Antibody**Catalog #: AMRe03266**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 41 kDa; Observed MW: 41 kDa

Antigen Information

Gene Name	HLA-A
Alternative Names	Aw-68; HLA class I histocompatibility antigen; A-28 alpha chain; MHC class I antigen A*68; HLA-A; MHC class I antigen HLA A heavy chain
Gene ID	3105.0
SwissProt ID	P04439
Immunogen	

Background


Major histocompatibility complex (MHC) molecules form an integral part of the immune response system. They are cell-surface

receptors that bind peptides and present them to T lymphocytes. HLA-A, -B and -C encode membrane anchored heavy chains which heterodimerize with a light chain (b-2-Microglobulin) to form MHC-I. Polymorphisms yield hundreds of HLA-A, -B and -C alleles.


Research Area

Immunology

Image Data



Western blot analysis of MHC class I in HeLa, A549 lysates using MHC Class I antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using MHC class I antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.