

**Product Name: Acetyl-Histone H4 (Lys5) Rabbit
Monoclonal Antibody
Catalog #: AMRe03268**

Summary

Production Name	Acetyl-Histone H4 (Lys5) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF,IP
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Acetylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	H4C1
Alternative Names	H4K5ac; H4F2; HIST; HIST1H4A; Histone H4.AC-H4K5; H4K5
Gene ID	121504
SwissProt ID	P62805.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200 IP: 1:20
Molecular Weight	Calculated MW: 11 kDa; Observed MW: 11 kDa

**Product Name: Acetyl-Histone H4 (Lys5) Rabbit
Monoclonal Antibody
Catalog #: AMRe03268**




Background

The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a member of the histone H4 family.


Research Area

Epigenetics and Nuclear Signaling

Image Data




Immunocytochemistry analysis of Acetyl-Histone H4 (Lys5) (green) in HeLa using Acetyl-Histone H4 (Lys5) antibody, and DAPI (blue).




Western blot analysis of Acetyl-Histone H4 (Lys5) in 3T3 lysates using Acetyl-Histone H4 (Lys5) antibody. Observed band size: 11 kDa.

**Product Name: Acetyl-Histone H4 (Lys5) Rabbit
Monoclonal Antibody
Catalog #: AMRe03268**



Western blot analysis of Acetyl-Histone H4 (Lys5) in A549, HL-60, U2OS, U87-MG lysates using Acetyl-Histone H4 (Lys5) antibody



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Acetyl-Histone H4 (Lys5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.