

**Product Name: Phospho-PAK4/5/6
(Ser474/Ser560/Ser602) Rabbit Monoclonal Antibody**
Catalog #: AMRe03279

Summary

Production Name	Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	PAK4/PAK5/PAK6
Alternative Names	PAK4; PAK5; PAK6
Gene ID	10298/57144/56924
SwissProt ID	O96013/Q9P286/Q9NQU5.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 64,80,75 kDa; Observed MW: 64,80,75 kDa

Background


**Product Name: Phospho-PAK4/5/6
(Ser474/Ser560/Ser602) Rabbit Monoclonal Antibody
Catalog #: AMRe03279**

Serine/threonine protein kinase that plays a role in a variety of different signaling pathways including cytoskeleton regulation, cell migration, growth, proliferation or cell survival.

Research Area

Signal Transduction

Image Data



Western blot analysis of Phospho-PAK4 + PAK5 + PAK6 (Ser474 + Ser560 + Ser602) in HeLa, A549, HL-60, U2OS, U87-MG lysates using Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) antibody.

Note

For research use only.