

Product Name: PI3 Kinase p85 alpha Rabbit Monoclonal Antibody
Catalog #: AMRe03743

Summary

Production Name	PI3 Kinase p85 alpha Rabbit Monoclonal Antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC/IF,IP
Reactivity	Human,Rat,Hamster

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	PIK3R1
Alternative Names	PIK3R1; GRB1; Phosphatidylinositol 3-kinase regulatory subunit alpha; PI3-kinase regulatory subunit alpha; PI3K regulatory subunit alpha; PtdIns-3-kinase regulatory subunit alpha; Phosphatidylinositol 3-kinase 85 kDa regulatory subunit alpha; PI3-kinase subunit p85-alpha; PtdIns-3-kinase regulatory subunit p85-alpha
Gene ID	5295
SwissProt ID	P27986.

Application

Dilution Ratio	WB 1:500-1:1000,ICC/IF 1:50-1:200,IP 1:20-1:50
-----------------------	--

Product Name: PI3 Kinase p85 alpha Rabbit Monoclonal Antibody
Catalog #: AMRe03743

Molecular Weight

Calculated MW: 84 kDa; Observed MW: 84 kDa


Background

Binds to activated (phosphorylated) protein-Tyr kinases, through its SH2 domain, and acts as an adapter, mediating the association of the p110 catalytic unit to the plasma membrane. Necessary for the insulin-stimulated increase in glucose uptake and glycogen synthesis in insulin-sensitive tissues. Plays an important role in signaling in response to FGFR1, FGFR2, FGFR3, FGFR4, KITLG/SCF, KIT, PDGFRA and PDGFRB. Likewise, plays a role in ITGB2 signaling.


Research Area

Signal Transduction

Image Data



Western blot analysis of PI3 Kinase p85 alpha in Jurkat, Rat Brain, C6, CHO-K1 lysates using PI3 Kinase p85 alpha antibody.



Immunocytochemistry analysis of PI3 Kinase p85 alpha (green) in MCF-7 using PI3 Kinase p85 alpha antibody, and DAPI (blue)

Note



**Product Name: PI3 Kinase p85 alpha Rabbit Monoclonal
Antibody
Catalog #: AMRe03743**

For research use only.