

Product Name: AKT Rabbit Monoclonal Antibody
Catalog #: AMRe03759



Summary

Production Name	AKT Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF,IP
Reactivity	Human,Mouse,Rat,Hamster

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	AKT1/AKT2/AKT3
Alternative Names	AKT1
Gene ID	207/208/10000
SwissProt ID	P31749/P31751/Q9Y243.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200 IP: 1:20
Molecular Weight	Calculated MW: 56 kDa; Observed MW: 56 kDa

Background

Product Name: AKT Rabbit Monoclonal Antibody
Catalog #: AMRe03759




Akt, also referred to as PKB or Rac, plays a critical role in controlling survival and apoptosis. This protein kinase is activated by insulin and various growth and survival factors to function in a wortmannin-sensitive pathway involving PI3 kinase. Akt is activated by phospholipid binding and activation loop phosphorylation at Thr308 by PDK1 and by phosphorylation within the carboxy terminus at Ser473.


Research Area

Signal Transduction


Image Data



Western blot analysis of AKT in 3T3, mouse lung lysates using AKT antibody.



Western blot analysis of AKT1/2/3 in Jurkat, rat Brain, C6, CHO-K1, Hela lysates using AKT antibody.



Product Name: AKT Rabbit Monoclonal Antibody
Catalog #: AMRe03759



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using AKT1/2/3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.