

Product Name: PFKFB3 Rabbit Monoclonal Antibody**Catalog #: AMRe03790**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.65mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 60 kDa; Observed MW: 60 kDa

Antigen Information

Gene Name	PFKFB3
Alternative Names	6-bisphosphatase; 6-P2ase 3; Fructose-2; IPFK2; PFK/FBPase 3; PFK2; PFKFB3; uPFK 2
Gene ID	5209
SwissProt ID	Q16875
Immunogen	A synthetic peptide corresponding to target protein


Background

Synthesis and degradation of fructose 2,6-bisphosphate.


Research Area

Cardiovascular

Image Data



Western blot analysis of PFKFB3 in Raw264.7, C6 lysates using PFKFB3 antibody.



Western blot analysis of PFKFB3 in K562 lysates using PFKFB3 antibody.