
Product Name: GAA Rabbit Monoclonal Antibody**Catalog #: AMRe03946**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW:105 kDa;Observed MW: 76 ,95,110kDa

Antigen Information

Gene Name	GAA
Alternative Names	Acid maltase; GAA; Glucosidase alpha acid; LYAG
Gene ID	2548
SwissProt ID	P10253
Immunogen	A synthetic peptide of human GAA


Background

Essential for the degradation of glycogen to glucose in lysosomes.

Research Area

Signal Transduction

Image Data



Western blot analysis of GAA in 293T, Human placenta lysates using GAA antibody.