

Product Name: CAD Rabbit Monoclonal Antibody**Catalog #: AMRe03951**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW:243 kDa;Observed MW: 243 kDa

Antigen Information

Gene Name	CAD
Alternative Names	CAD; CAD protein
Gene ID	790
SwissProt ID	P27708
Immunogen	A synthetic peptide of human CAD

Background

Carbamoyl phosphate synthetase-aspartate carbamoyltransferase-dihydroorotase (CAD) is a multifunctional protein that initiates and regulates mammalian de novo pyrimidine biosynthesis.

Research Area

Cell Biology

Image Data

Western blot analysis of CAD in 293 lysates using CAD antibody.

