

Product Name: Yes1 Rabbit Monoclonal Antibody**Catalog #: AMRe03958**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW:61 kDa;Observed MW: 61 kDa

Antigen Information

Gene Name	YES1
Alternative Names	YES1; YES; Tyrosine-protein kinase Yes; Proto-oncogene c-Yes; p61-Yes
Gene ID	7525
SwissProt ID	P07947
Immunogen	A synthetic peptide of human Yes1

Background

Promotes infectivity of *Neisseria gonorrhoeae* in epithelial cells by phosphorylating MCP/CD46.

Research Area

Signal Transduction

Image Data

Western blot analysis of Yes1 in SW480, Jurkat lysates using Yes1 antibody.

