

Product Name: C1orf93 Rabbit Monoclonal Antibody**Catalog #: AMRe03963**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1.21mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW:21 kDa;Observed MW: 21 kDa

Antigen Information

Gene Name	PRXL2B
Alternative Names	C1orf93
Gene ID	127281
SwissProt ID	Q8TBF2
Immunogen	Recombinant protein of human C1orf93

Background

Catalyzes the reduction of prostaglandin-ethanolamide H2 (prostamide H2) to prostamide F(2 α) with NADPH as proton donor.

Research Area

Cell Biology

Image Data

Western blot analysis of C1orf93 in MCF7 lysates using C1orf93 antibody.

