
Product Name: Syntaxin 3 Rabbit Monoclonal Antibody**Catalog #: AMRe03964**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
Molecular Weight	Calculated MW:33 kDa;Observed MW: 33 kDa

Antigen Information

Gene Name	STX3
Alternative Names	Stx3; STX3A; Syntaxin 3A
Gene ID	6809
SwissProt ID	Q13277
Immunogen	A synthetic peptide of human Syntaxin 3


Background

Potentially involved in docking of synaptic vesicles at presynaptic active zones.

Research Area

Neuroscience

Image Data



Western blot analysis of Syntaxin 3 in HepG2, Raw 264.7 lysates using Syntaxin 3 antibody.