

Product Name: INTS10 Rabbit Monoclonal Antibody**Catalog #: AMRe03965**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW:82 kDa;Observed MW: 82 kDa

Antigen Information


Gene Name	INTS10
Alternative Names	integrator complex subunit 10; INT10; C8orf35
Gene ID	55174
SwissProt ID	Q9NVR2
Immunogen	A synthetic peptide of human INTS10

Background

Component of the Integrator (INT) complex, a complex involved in the small nuclear RNAs (snRNA) U1 and U2 transcription and in their 3'-box-dependent processing.

Research Area

Image Data



Western blot analysis of INTS10 in HeLa, Jurkat lysates using INTS10 antibody.