
Product Name: PI 3 Kinase p150 Rabbit Monoclonal Antibody**Catalog #: AMRe03969**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC/IF
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,ICC/IF 1:50-1:200
Molecular Weight	Calculated MW:153 kDa;Observed MW: 153 kDa

Antigen Information

Gene Name	PIK3R4
Alternative Names	P150; Phosphoinositide 3 kinase regulatory subunit 4; PI3 kinase p150 subunit; PIK3R4; VPS15
Gene ID	30849
SwissProt ID	Q99570
Immunogen	A synthetic peptide of human PIK3R4


Background

Regulatory subunit of the PI3K complex. May regulate membrane trafficking late in the endocytic pathway.

Research Area

Signal Transduction

Image Data



Western blot analysis of PI 3 Kinase p150 in Jurkat,Raji lysates using PI 3 Kinase p150 antibody.