

Summary

Production Name	Granulin Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Purification

Immunogen

Gene Name	GRN
Alternative Names	Acrogranin; CLN11; GEP; GP88; Granulins; GRN; PCDGF; PEPI; PGRN; Proepithelin; Progranulin
Gene ID	2896.0
SwissProt ID	P28799.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW:64 kDa;Observed MW: 55,88 kDa

Background

Product Name: Granulin Rabbit Monoclonal Antibody
Catalog #: AMRe03970




Granulins have possible cytokine-like activity. They may play a role in inflammation, wound repair, and tissue remodeling.

Research Area

Neuroscience

Image Data



Western blot analysis of Granulin in 293T lysates using Granulin antibody.

Note

For research use only.