

Product Name: GARNL1 Rabbit Monoclonal Antibody
Catalog #: AMRe03975



Summary

Production Name	GARNL1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Purification

Immunogen

Gene Name	RALGAPA1
Alternative Names	p240; GRIPE; GARNL1; TULIP1; RaGAPalpha1
Gene ID	253959
SwissProt ID	Q6GYQ0.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200
Molecular Weight	Calculated MW:230 kDa;Observed MW: 230 kDa

Background

Catalytic subunit of the heterodimeric RaGAP1 complex which acts as a GTPase activator for the Ras-like small GTPases

Product Name: GARNL1 Rabbit Monoclonal Antibody
Catalog #: AMRe03975




RALA and RALB.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of GARNL1 in 293, lysates using GARNL1 antibody.

Note

For research use only.