

Product Name: HOXA9 Rabbit Monoclonal Antibody**Catalog #: AMRe03985**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW:30 kDa;Observed MW: 30 kDa

Antigen Information

Gene Name	HOXA9
Alternative Names	HOX1; ABD-B; HOX1G; HOX1.7
Gene ID	3205
SwissProt ID	P31269
Immunogen	A synthetic peptide of human HOXA9

Background

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Research Area

Epigenetics and Nuclear Signaling

Image Data

Western blot analysis of HOXA9 in U937 lysates using HOXA9 antibody.

