
Product Name: Lyn Rabbit Monoclonal Antibody**Catalog #: AMRe03992**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW:59 kDa;Observed MW: 59 kDa

Antigen Information

Gene Name	LYN
Alternative Names	LYN; JTK8; Tyrosine-protein kinase Lyn; Lck/Yes-related novel protein tyrosine kinase; V-yes-1 Yamaguchi sarcoma viral related oncogene homolog; p53Lyn; p56Lyn
Gene ID	4067
SwissProt ID	P07948
Immunogen	A synthetic peptide corresponding to target protein

Background

Lyn a tyrosine kinase of the Src family that plays a role in B cell antigen-receptor signaling.

Research Area

Signal Transduction

Image Data

Western blot analysis of Lyn in HeLa lysates using Lyn antibody.

