

Product Name: CKS2 Rabbit Monoclonal Antibody**Catalog #: AMRe03996**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.3mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW:10 kDa;Observed MW: 10 kDa

Antigen Information

Gene Name	CKS2
Alternative Names	CKS-2,Cyclin-dependent kinases regulatory subunit 2
Gene ID	1164.0
SwissProt ID	P33552
Immunogen	A synthetic peptide of human CKS2

Background

Binds to the catalytic subunit of the cyclin dependent kinases and is essential for their biological function.

Research Area

Cell Biology

Image Data

Western blot analysis of CKS2 in Caco-2, 293T lysates using CKS2 antibody.

