

Product Name: FGL1 Rabbit Monoclonal Antibody**Catalog #: AMRe04009**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW:36 kDa;Observed MW: 36 kDa

Antigen Information

Gene Name	FGL1
Alternative Names	FGL1; HFREP1; Fibrinogen-like protein 1; HP-041; Hepassocin; Hepatocyte-derived fibrinogen-related protein 1; HFREP-1; Liver fibrinogen-related protein 1; LFIRE-1
Gene ID	2267
SwissProt ID	Q08830
Immunogen	A synthetic peptide of human FGL1

Background

Immune suppressive molecule that inhibits antigen-specific T-cell activation by acting as a major ligand of LAG3.

Research Area

Cardiovascular

Image Data

Western blot analysis of FGL1 in Human fetal liver lysates using FGL1 antibody.

