

Product Name: CD239/BCAM Rabbit Monoclonal Antibody**Catalog #: AMRe04011**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1.36mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
Molecular Weight	Calculated MW:67 kDa;Observed MW: 85kDa

Antigen Information

Gene Name	BCAM BCAM; LU; MSK19; Basal cell adhesion molecule; Auberger B antigen; B-CAM cell surface
Alternative Names	glycoprotein; F8/G253 antigen; Lutheran antigen; Lutheran blood group glycoprotein; CD239
Gene ID	4059
SwissProt ID	P50895
Immunogen	A synthetic peptide of human CD239/BCAM

Background

Laminin alpha-5 receptor. May mediate intracellular signaling.

Research Area

Cardiovascular

Image Data

Western blot analysis of CD239/BCAM in HeLa lysates using CD239/BCAM antibody.

