

Product Name: Plexin A1 Rabbit Monoclonal Antibody**Catalog #: AMRe04017**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,FC
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,FC 1:50-1:100
Molecular Weight	Calculated MW:211 kDa;Observed MW: 211 kDa

Antigen Information

Gene Name	PLXNA1
Alternative Names	NOV; NOV; Plexin-A1; Plexin1; PlexinA1; PLXN1; Plxna1
Gene ID	5361
SwissProt ID	Q9UIW2
Immunogen	A synthetic peptide of human Plexin A1

Background

Plays a role in axon guidance, invasive growth and cell migration.

Research Area

Neuroscience

Image Data

Western blot analysis of Plexin A1 in K562 lysates using Plexin A1 antibody.

