

Product Name: Dishevelled 3 Rabbit Monoclonal Antibody**Catalog #: AMRe04028**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,FC
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100
Molecular Weight	Calculated MW:78 kDa;Observed MW: 95 kDa

Antigen Information

Gene Name	DVL3
Alternative Names	DRS3,KIAA0208,Dishevelled-3,DSH homolog 3,Segment polarity protein dishevelled homolog DVL-3
Gene ID	1857.0
SwissProt ID	Q92997
Immunogen	A synthetic peptide of human Dishevelled 3

Background

Involved in the signal transduction pathway mediated by multiple Wnt genes.

Research Area

Neuroscience

Image Data

Western blot analysis of Dishevelled 3 in K562 lysates using Dishevelled 3 antibody.

