

Product Name: FGF2 Rabbit Monoclonal Antibody
Catalog #: AMRe04037



Summary

Production Name	FGF2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	FGF2
Alternative Names	FGF2; FGFB; Fibroblast growth factor 2; FGF-2; Basic fibroblast growth factor; bFGF; Heparin-binding growth factor 2; HBGF-2
Gene ID	2247
SwissProt ID	P09038.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 31 kDa; Observed MW: 19 kDa

Background

Product Name: FGF2 Rabbit Monoclonal Antibody
Catalog #: AMRe04037




Plays an important role in the regulation of cell survival, cell division, angiogenesis, cell differentiation and cell migration.
Functions as potent mitogen in vitro. Can induce angiogenesis.

Research Area

Cardiovascular

Image Data



Western blot analysis of FGF2 in HeLa lysates using FGF2 antibody.

Note

For research use only.