

Product Name: DAZAP1 Rabbit Monoclonal antibody**Catalog #: AMRe04071**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 43 kDa; Observed MW: 43 kDa

Antigen Information

Gene Name	DAZAP1
Alternative Names	Deleted in azoospermia-associated protein 1
Gene ID	26528
SwissProt ID	Q96EP5
Immunogen	Recombinant protein of human DAZAP1


Background

DAZ-associated protein 1

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of DAZAP1 in K562, C6, 3T3, HeLa lysates using DAZAP1 antibody.