

Product Name: eIF3B Rabbit Monoclonal antibody
Catalog #: AMRe04078



Summary

Production Name	eIF3B Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF, IP
Reactivity	Human, Mouse, Monkey

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	EIF3B
Alternative Names	eIF3b; EIF3S9; hPrt1; PRT1
Gene ID	8662
SwissProt ID	P55884.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200 IP: 1:20
Molecular Weight	Calculated MW: 92 kDa; Observed MW: 117 kDa

Background

Product Name: eIF3B Rabbit Monoclonal antibody
Catalog #: AMRe04078




Component of the eukaryotic translation initiation factor 3 (eIF-3) complex, which is required for several steps in the initiation of protein synthesis.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of eIF3B in Raji, COS7, 3T3, HeLa lysates using eIF3B antibody.



Immunocytochemistry analysis of eIF3B(green) in HeLa using eIF3B antibody, and DAPI(blue)

Note

For research use only.