

Product Name: Flotillin 1 Rabbit Monoclonal Antibody
Catalog #: AMRe04079



Summary

Production Name	Flotillin 1 Rabbit Monoclonal Antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF, IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	FLOT1
Alternative Names	Flotillin-1
Gene ID	10211
SwissProt ID	O75955.

Application

Dilution Ratio	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
Molecular Weight	Calculated MW: 47 kDa; Observed MW: 47 kDa

Background

Product Name: Flotillin 1 Rabbit Monoclonal Antibody
Catalog #: AMRe04079




May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles.


Research Area

Cell Biology

Image Data



Western blot analysis of Flotillin 1 in K562, HeLa lysates using Flotillin 1 antibody.



Immunocytochemistry analysis of Flotillin 1(green) in HeLa using Flotillin 1 antibody, and DAPI(blue)

Note

For research use only.