

Product Name: ICAD Rabbit Monoclonal Antibody
Catalog #: AMRe04085



Summary

Production Name	ICAD Rabbit Monoclonal Antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	DFFA
Alternative Names	DFFA; DFF1; DFF45; H13; DNA fragmentation factor subunit alpha; DNA fragmentation factor 45 kDa subunit; DFF-45; Inhibitor of CAD; ICAD
Gene ID	1676
SwissProt ID	O00273.

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 37 kDa; Observed MW: 45 kDa

Product Name: ICAD Rabbit Monoclonal Antibody
Catalog #: AMRe04085




Background

Inhibitor of the caspase-activated DNase (DFF40).

Research Area

Cell Biology

Image Data



Western blot analysis of ICAD in rat Brain, MCF-7, C6 lysates using ICAD antibody.

Note

For research use only.