

Product Name: MCM7 Rabbit Monoclonal antibody
Catalog #: AMRe04090



Summary

Production Name	MCM7 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	MCM7
Alternative Names	MCM7; CDC47; MCM2; DNA replication licensing factor MCM7; CDC47 homolog; P1.1-MCM3
Gene ID	4176
SwissProt ID	P33993.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 81 kDa; Observed MW: 81 kDa

Background

Product Name: MCM7 Rabbit Monoclonal antibody
Catalog #: AMRe04090




Acts as component of the MCM2-7 complex (MCM complex) which is the putative replicative helicase essential for 'once per cell cycle' DNA replication initiation and elongation in eukaryotic cells. The active ATPase sites in the MCM2-7 ring are formed through the interaction surfaces of two neighboring subunits such that a critical structure of a conserved arginine finger motif is provided in trans relative to the ATP-binding site of the Walker A box of the adjacent subunit.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of MCM7 in K562, C6, 3T3, HeLa lysates using MCM7 antibody.

Note

For research use only.