

Product Name: MLKL Rabbit Monoclonal antibody
Catalog #: AMRe04091



Summary

Production Name	MLKL Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	MLKL
Alternative Names	MLKL
Gene ID	197259
SwissProt ID	Q8NB16.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 54 kDa; Observed MW: 54 kDa

Background

Product Name: MLKL Rabbit Monoclonal antibody
Catalog #: AMRe04091




Pseudokinase that plays a key role in TNF-induced necroptosis, a programmed cell death process. Activated following phosphorylation by RIPK3, leading to homotrimerization, localization to the plasma membrane and execution of programmed necrosis characterized by calcium influx and plasma membrane damage. Does not have protein kinase activity.


Research Area

Signal Transduction

Image Data



Western blot analysis of MLKL in HT-29 lysates using MLKL antibody.



Western blot analysis of MLKL in U251MG, Hela lysates using MLKL antibody.

Note

For research use only.