

Summary

Production Name	PRMT1 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	PRMT1
Alternative Names	ANM 1; ANM1; ANM1_HUMAN; complete cds; HCP 1; HCP1; Heterogeneous nuclear ribonucleoprotein methyltransferase 1 like 2; Heterogeneous nuclear ribonucleoproteins methyltransferase like 2; Heterogeneous nuclear ribonucleoproteins methyltransferase like2; Histone-arginine N-methyltransferase PRMT1; HMT 2; HMT1 (hnRNP methyltransferase; HMT1 (hnRNP methyltransferase S. cerevisiae) like 2; HMT1 hnRNP methyltransferase; HMT1 hnRNP methyltransferase like 2 (S. cerevisiae); HMT1 hnRNP methyltransferase like 2; HMT1 hnRNP methyltransferase-like 2 (S. cerevisiae); HMT2; HRMT1 L2; HRMT1L 2; HRMT1L2; Human mRNA for suppressor for yeast mutant; Human mRNA for suppressor for yeast mutant complete cds; Interferon receptor 1 bound protein 4; Interferon receptor 1 bound protein4; Interferon receptor 1-bound protein 4; Interferon receptor 1bound

Product Name: PRMT1 Rabbit Monoclonal antibody
Catalog #: AMRe04103



	protein 4; IR1 B4; IR1B 4; IR1B4; Mrmt 1; Mrmt1; PRMT 1; PRMT1; Protein arginine methyltransferase 1; Protein arginine N methyltransferase 1; Protein arginine N methyltransferase1; Protein arginine N-methyltransferase 1; R1B4; S. cerevisiae like 2.
Gene ID	3276
SwissProt ID	Q99873.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW: 42 kDa; Observed MW: 42 kDa


Background

Arginine methyltransferase that methylates (mono and asymmetric dimethylation) the guanidino nitrogens of arginyl residues present in proteins such as ESR1, histone H2, H3 and H4, PIAS1, HNRNPA1, HNRNPD, NFATC2IP, SUPT5H, TAF15 and EWS. Constitutes the main enzyme that mediates monomethylation and asymmetric dimethylation of histone H4 'Arg-4' (H4R3me1 and H4R3me2a, respectively), a specific tag for epigenetic transcriptional activation.

Research Area


Cell Biology

Image Data



Western blot analysis of PRMT1 in K562, C6 lysates using PRMT1 antibody.

Product Name: PRMT1 Rabbit Monoclonal antibody
Catalog #: AMRe04103



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using PRMT1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.