

Product Name: Integrin alpha 5 Rabbit Monoclonal antibody
Catalog #: AMRe04142

Summary

Production Name	Integrin alpha 5 Rabbit Monoclonal antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	ITGA5
Alternative Names	ITGA5; FNRA; Integrin alpha-5; CD49 antigen-like family member E; Fibronectin receptor subunit alpha; Integrin alpha-F; VLA-5; CD49e
Gene ID	3678
SwissProt ID	P08648.

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
Molecular Weight	Calculated MW: 115 kDa; Observed MW: 150 kDa

Product Name: Integrin alpha 5 Rabbit Monoclonal antibody
Catalog #: AMRe04142




Background

Integrin alpha-5/beta-1 is a receptor for fibronectin and fibrinogen. It recognizes the sequence R-G-D in its ligands. ITGA5:ITGB1 binds to PLA2G2A via a site (site 2) which is distinct from the classical ligand-binding site (site 1) and this induces integrin conformational changes and enhanced ligand binding to site 1. ITGA5:ITGB1 acts as a receptor for fibrillin-1 (FBN1) and mediates R-G-D-dependent cell adhesion to FBN1.


Research Area

Signal Transduction

Image Data



Western blot analysis of Integrin alpha 5 in K562, 3T3 lysates using Integrin alpha 5 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Integrin alpha 5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.