

**Product Name: Phospho-IRF3 (Ser386) (9P16) Rabbit
Monoclonal Antibody
Catalog #: AMRe05928**



Summary

Production Name	Phospho-IRF3 (Ser386) (9P16) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
Purification	Affinity purification

Immunogen

Gene Name	IRF3
Alternative Names	
Gene ID	3661.0
SwissProt ID	Q14653.

Application

Dilution Ratio	WB 1:1000-1:5000
Molecular Weight	47 kDa


Background

**Product Name: Phospho-IRF3 (Ser386) (9P16) Rabbit
Monoclonal Antibody
Catalog #: AMRe05928**



Research Area

Image Data



Western blot analysis of Phospho-IRF3 (S386) in MCF-7 lysates using Phospho-IRF3 (Ser386) antibody.

Note

For research use only.