

**Product Name: Phospho-POLR2A (S2) (1B7) Rabbit
Monoclonal Antibody
Catalog #: AMRe05979**

Summary

Production Name	Phospho-POLR2A (S2) (1B7) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,FC,IP,IF-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Immunogen

Gene Name	POLR2A
Alternative Names	POLR2A; POLR2; RNA polymerase II CTD repeat YSPTSPS;
Gene ID	5430.0
SwissProt ID	P24928. A synthetic phosphopeptide corresponding to residues surrounding Ser2 of human RNA polymerase II CTD repeat YSPTSPS

Application


Dilution Ratio	WB 1:1000, IHC-P/IF-P 1:200-1:2000, ICC/IF 1:200-1:500, FCM 1:200-1:500, IP 1:20-1:50
Molecular Weight	192kDa

Background

During transcription elongation, Pol II moves on the template as the transcript elongates. Elongation is influenced by the phosphorylation status of the C-terminal domain (CTD) of Pol II largest subunit (RPB1), which serves as a platform for assembly of factors that regulate transcription initiation, elongation, termination and mRNA processing. DNA-dependent RNA polymerase catalyzes the transcription of DNA into RNA using the four ribonucleoside triphosphates as substrates. Largest and catalytic component of RNA polymerase II which synthesizes mRNA precursors and many functional non-coding RNAs. Forms the polymerase active center together with the second largest subunit. Pol II is the central component of the basal RNA polymerase II transcription machinery. It is composed of mobile elements that move relative to each other. RPB1 is part of the core element with the central large cleft, the clamp element that moves to open and close the cleft and the jaws that are thought to grab the incoming DNA template. At the start of transcription, a single-stranded DNA template strand of the promoter is positioned within the central active site cleft of Pol II. A bridging helix emanates from RPB1 and crosses the cleft near the catalytic site and is thought to promote translocation of Pol II by acting as a ratchet that moves the RNA-DNA hybrid through the active site by switching from straight to bent conformations at each step of nucleotide addition. During transcription elongation, Pol II moves on the template as the transcript elongates. Elongation is influenced by the phosphorylation status of the C-terminal domain (CTD) of Pol II largest subunit (RPB1), which serves as a platform for assembly of factors that regulate transcription initiation, elongation, termination and mRNA processing. Regulation of gene expression levels depends on the balance between methylation and acetylation levels of the CTD-lysines (By similarity). Initiation or early elongation steps of transcription of growth-factors-induced immediate early genes are regulated by the acetylation status of the CTD (PubMed: <http://www.uniprot.org/citations/24207025> target="_blank">24207025). Methylation and dimethylation have a repressive effect on target genes expression (By similarity).

Research Area

Image Data



**Product Name: Phospho-POLR2A (S2) (1B7) Rabbit
Monoclonal Antibody
Catalog #: AMRe05979**



Western blot analysis of extracts from HeLa cells using RM5247 at 1:1000.

Note

For research use only.