

Product Name: Phospho-Smad3 (Ser423/425) (12O16) Rabbit Monoclonal Antibody
Catalog #: AMRe06012

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:100-1:200,ICC/IF 1:100-1:200
Molecular Weight	48kDa

Antigen Information


Gene Name	SMAD3
Alternative Names	
Gene ID	4088.0
SwissProt ID	P84022
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser423/Ser425 of human Smad3

Background

Research Area

Signal Transduction

Image Data



Western blot detection of Phospho-Smad3 (Ser423/425) in Jurkat,C6,CHO-K1,HeLa cell lysates using Phospho-Smad3 (Ser423/425) antibody(1:1000 diluted).