

Product Name: AGTR2 (2I6) Rabbit Monoclonal Antibody
Catalog #: AMRe06686

Summary

Production Name	AGTR2 (2I6) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type
Buffer	preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	AGTR2
Alternative Names	AGTR2; Angiotensin ii receptor type 2; Angiotensin II type-2 receptor; AT2 receptor; AT2; At2r; MRX88; Angiotensin ii at2; Angiotensin receptor 2; Type-2 angiotensin II receptor; ATGR2;
Gene ID	186.0
SwissProt ID	P50052.

Application

Dilution Ratio	WB 1:1000-1:5000, IP 1:10-1:100
Molecular Weight	41kDa


Product Name: AGTR2 (2I6) Rabbit Monoclonal Antibody
Catalog #: AMRe06686

Background

Receptor for angiotensin II. Cooperates with MTUS1 to inhibit ERK2 activation and cell proliferation. Receptor for angiotensin II. Cooperates with MTUS1 to inhibit ERK2 activation and cell proliferation.

Research Area

Image Data



Western blot analysis of AGTR2 expression in HepG2 cell lysate.

Note

For research use only.