
Product Name: AMPK alpha (14Y14) Rabbit Monoclonal Antibody**Catalog #: AMRe06840**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IP 1:20-1:50
Molecular Weight	64,62kDa

Antigen Information


Gene Name	PRKAA2
Alternative Names	
Gene ID	5563.0
SwissProt ID	P54646
Immunogen	Recombinant protein of human AMPK alpha 2

Background

Research Area

Neuroscience

Image Data



Western blot detection of AMPK alpha 1/2 in Jurkat cell lysates using AMPK alpha 1/2 antibody(1:500 diluted).