

Product Name: Annexin A10 (14O15) Rabbit Monoclonal Antibody
Catalog #: AMRe06915

Summary

Production Name	Annexin A10 (14O15) Rabbit Monoclonal Antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IF-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	ANXA10
Alternative Names	annexin 14; annexin A10; annexin-14; ANX10; ANX14;
Gene ID	11199.0
SwissProt ID	Q9UJ72.

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:200-1:2000,IF-P 1:200-1:2000
Molecular Weight	37kDa


Product Name: Annexin A10 (14O15) Rabbit Monoclonal Antibody
Catalog #: AMRe06915

Background

Belongs to the annexin family. Contains 4 annexin repeats.

Research Area

Image Data



Western blot analysis of Annexin A10 expression in HeLa cell lysate.

Note

For research use only.