

Product Name: Aquaporin 2 (13E19) Rabbit Monoclonal Antibody**Catalog #: AMRe07079**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP,IF-P
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:500-1:2000,IP 1:20-1:50,IF-P 1:500-1:2000
Molecular Weight	29kDa

Antigen Information


Gene Name	AQP2
Alternative Names	
Gene ID	359.0
SwissProt ID	P41181
Immunogen	A synthetic peptide of human Aquaporin 2

Background

Research Area

Cell Biology

Image Data



Western blot detection of Aquaporin 2 in Raoms,Raji cell lysates using Aquaporin 2 antibody(1:1000 diluted).