

---

**Product Name: ATOH1 (17T16) Rabbit Monoclonal Antibody****Catalog #: AMRe07313**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	0.5mg/ml. The concentration of this product may be batch-dependent.
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
<b>Purification</b>	Affinity purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:2000
<b>Molecular Weight</b>	38kDa

**Antigen Information**

<b>Gene Name</b>	ATOH1
<b>Alternative Names</b>	ATH1; ATOH1; bHLHa14; hATH1; MATH1;
<b>Gene ID</b>	474.0
<b>SwissProt ID</b>	Q92858
<b>Immunogen</b>	A synthetic peptide of human MATH1/HATH1


**Background**

Activates E box-dependent transcription in collaboration with TCF3/E47, but the activity is completely antagonized by the

negative regulator of neurogenesis HES1. May play a role in the differentiation of subsets of neural cells by activating E box-dependent transcription. Transcriptional regulator. Activates E box-dependent transcription in collaboration with TCF3/E47, but the activity is completely antagonized by the negative regulator of neurogenesis HES1. Plays a role in the differentiation of subsets of neural cells by activating E box-dependent transcription (By similarity).

## Research Area

## Image Data



Western blot analysis of MATH1/HATH1 expression in HepG2 cell lysate.