

Product Name: Biglycan (9M18) Rabbit Monoclonal Antibody
Catalog #: AMRe07560



Summary

Production Name	Biglycan (9M18) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,FC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type
Buffer	preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	BGN
Alternative Names	BGN; Biglycan; DSPG1; PG S1; PGI; SEMDX; SLRR1A;
Gene ID	633.0
SwissProt ID	P21810.

Application

Dilution Ratio	WB 1:1000-1:5000, FCM 1:100-1:500
Molecular Weight	42kDa

Product Name: Biglycan (9M18) Rabbit Monoclonal Antibody
Catalog #: AMRe07560




Background

May be involved in collagen fiber assembly. May be involved in collagen fiber assembly.

Research Area

Image Data



Western blot analysis of Biglycan expression in C2C12 cell lysate.

Note

For research use only.