

Product Name: Calponin (6N15) Rabbit Monoclonal Antibody
Catalog #: AMRe07871



Summary

Production Name	Calponin (6N15) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,IF-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type
Buffer	preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	CNN1
Alternative Names	Calponin-1; Basic calponin; Calponin H1, smooth muscle; CNN1
Gene ID	1264.0
SwissProt ID	P51911.

Application

Dilution Ratio	WB 1:1000-1:5000, IHC-P/IF-P 1:200-1:1000, ICC/IF 1:200-1:500
Molecular Weight	33kDa

Product Name: Calponin (6N15) Rabbit Monoclonal Antibody
Catalog #: AMRe07871




Background

Calponin regulates smooth muscle cell contraction and is a marker of smooth muscle cell differentiation. Calponin, an Actin- and Tropomyosin-binding protein, is characterized as an inhibitory factor of smooth-muscle actomyosin activity. Calponin is implicated in the regulation of smooth muscle contraction through its interaction with F-Actin and inhibition of the Actin-activated MgATPase activity of phosphorylated Myosin. Thin filament-associated protein that is implicated in the regulation and modulation of smooth muscle contraction. It is capable of binding to actin, calmodulin and tropomyosin. The interaction of calponin with actin inhibits the actomyosin Mg-ATPase activity (By similarity).

Research Area

Image Data



Western blot analysis of Calponin expression in Mouse stomach lysate.

Note

For research use only.