

Product Name: CCR8 (11B19) Rabbit Monoclonal Antibody
Catalog #: AMRe08161

Summary

Production Name	CCR8 (11B19) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ICC/IF,FC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type
Buffer	preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	CCR8
Alternative Names	CY6; TER1; CCR-8; CKRL1; CDw198; CMKBR8; GPRCY6; CMKBRL2; CC-CKR-8; TER1;
Gene ID	1237.0
SwissProt ID	P51685.

Application

Dilution Ratio	WB 1:1000, ICC/IF 1:100, FCM 1:500-1:2000
Molecular Weight	41kDa


Product Name: CCR8 (11B19) Rabbit Monoclonal Antibody
Catalog #: AMRe08161

Background

Receptor for the chemokine CCL1/SCYA1/I-309. May regulate monocyte chemotaxis and thymic cell line apoptosis. Alternative coreceptor with CD4 for HIV-1 infection. Receptor for the chemokine CCL1/SCYA1/I-309. May regulate monocyte chemotaxis and thymic cell line apoptosis. Alternative coreceptor with CD4 for HIV-1 infection.

Research Area

Image Data



Western blot analysis of CCR8 expression in Jurkat cell lysate.

Note

For research use only.