
Product Name: CD226 (7Z13) Rabbit Monoclonal Antibody**Catalog #: AMRe08289**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,IP
Reactivity	Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:200-1:2000,ICC/IF 1:100-1:200,IP 1:20-1:50
Molecular Weight	38kDa

Antigen Information


Gene Name	Cd226
Alternative Names	
Gene ID	225825.0
SwissProt ID	Q8K4F0
Immunogen	Recombinant protein of mouse CD226

Background

Research Area

Immunology

Image Data



Western blot detection of CD226 in mouse spleen cell lysates using CD226 antibody(1:1000 diluted).